

List of commands for the USB interface

Parameters for the communication port :

Baud rate : 115200 Bauds Data bits : 8 Stop bits : 1
No parity Flow control : none Character for transmission end : LF (0xA)

List of the commands :

Important note: the commands have always 5 characters and it is important to respect the syntax of each of those (i.e. capital letters have to be used).

«HELP_» : (code version 5.7 and further)

This command allows to know all the information related to the different commands.

«_RAZ_» :

This commands sets all the outputs of the DAC to 0V.

« INFOS » :

This commands allows to know all the basic informations of the nanopositionning system (name of the stage, number of axis, maximum stroke, etc...)

« REOFF » : (code version 6.2 and further)

This function allows to disable the answers of the controller in order to reduce the latency. Warning: with this function the controller will not give answers until you reactivate it with RE_ON. The only command which will give an answer will be “INFOS” and “SERIE”.

« RE_ON » : (code version 6.2 and further)

This function allows you to reactivate the answers of the controller when it has been previously deactivated.

« **MOVEX** » :

This command allows to move the stage to an exact position. The same function MOVEY and MOVEZ exist for the Y axis and the Z axis.

For example, if one wants to move to the 200 μm position (ie 200000 nm), one must write :
« MOVEX 200u » or « MOVEX 200000n »

« **MOXYZ** » :

This function works in the same way as MOVEX but allows to make a movement on all axis at the same time.

ex: MOXYZ 100u 50u 25u

« **MOVRX** » : (code version 6.1 and further)

This command allows to move the X axis relatively to the current position. The same function MOVRY and MOVRZ exist for the Y axis and the Z axis. For example, if one wants to move 30 μm in the positive direction, one must write : « MOVRX +30u » ,
If one wants to move in the other direction, one must write « MOVRX -30u ».

« **MRXYZ** » : (code version 6.1 and further)

This function works in the same way as MOVRX but allows to realize a relative movement on all 3 axis at the same time.

ex: MRXYZ 10u 10u 10u

« **GET_X** » : This command allows to read the position of the X axis (similar function exist for the Y axis and the Z axis : GET_Y and GET_Z)

« **GEXYZ** » : (code version 5.7 and further)

This function works in the same way as GET_X but allows to acquire the position of the

3 axis at the same time.

« **STIME** » :

This command allows the modification of the time between 2 points which are sent by the RAM memory of the USB interface to the nanopositioner.

Consequently,

If one want to have a time between points of 200ms, one must write : « STIME 200m ».

If one want to have a time between points of 20µs, one must write « STIME 20u ».

If one want to have a time between points of 2s, one must write « STIME 2s ».

« **SHTIM** » :

This function allows you to set the shooting time of your acquisition when using the CHAIO function with the prefix t. The TTL pulses will then be set to the optimal position.

Consequently,

If I want to have a shooting time of 200ms, I must write : « SHTIM 200m ».

If I want to have a shooting time of 200µs, I must write « SHTIM 200u ».

If I want to have a shooting time of 1s, I must write « SHTIME 1s ».

« **GTIME** » :

This command allows you to know the time between each point.

« **SWF_X** »

This command allows you to set up a ramp waveform for the X axis. Example:

"SWF_X 100 0u 100u" loads a waveform of 100 equal sized steps from 0u to 100u.

The same function SWF_Y and SWF_Z exist for the Y axis and the Z axis.

« ARBWF »

This command prepares the controller for receiving arbitrary waveform data.

Data storage is allocated for the data when ARBWF is executed, data points must then be loaded using ADDPn.

Example:

« ARBWF 100 20 10 » This allocates storage for 100 X axis data points, 20 Y axis data points and 10 Z axis data points.

« ADDPX »

This command allows adding of data points to a waveform.

This command is only valid when used after ARBWF has been used to allocate data storage for the appropriate waveform data.

Example:

«ADDPX 100u» This command adds a point to the end of the X waveform data set. The command accepts the same format of position as MOVEX.

The same function ADDPY and ADDPZ exist for the Y axis and the Z axis.

« RUNWF »

This command allows to launch the scan defined by the function SWF_X, SWF_Y and SWF_Z. X axis is the first axis. Y axis is the second axis. Z axis is the third axis. One should write « RUNWF ».

« RUXYZ » (code version 5.7 and further)

This command is the same as RUNWF but also works for all axis in every sequences

Example: RUZXY, RUYXZ ,...

« RUXY_ » (code version 5.7 and further)

This command is the same as RUNWF but works for 2 axes in every sequences

Example: RUZX_, RUYX_ ,...

« RUX__» (code version 5.7 and further)

This command is the same as RUNWF but works only for 1 axe in all orders

Example: RUZ__, RUY__ ,...

« REXYZ» « REXY_» (code version 5.7 and further)

This command is the same as RUXYZ and RUXZ_ but with the first axis going back and forth.

« SWF_A » (code version 5.7 and further)

It is possible to run the RUXYZ several but defining an arbitrary value called A. The value you will give to A will allow to define the number of repetition of a RUXYZ function, by using the RXYZA function .

"RXYZA"(code version 5.7 and further)

This command allows to repeat several times

Example:

Example:

« SWF_A 10»” You set the number of repetition at 10.

« RXYZA » The function RUXYZ will be repeated 10 times.

«PAUSE» «PLAYY» and «STOPP» (code version 5.7 and further)

During a move with RUNWF, RUX__, RUXY_, RUXYZ or RUN3D, this function allows you to stop the function or pause it and start it again.

« ARB3D »

This command prepares the controller for receiving a sequence of arbitrary 3D locations. Data storage is allocated for the data when ARB3D is executed, data points must then be

loaded using ADD3D. The maximum number of 3D coordinates that can be allocated is about 2000.

Example: «ARB3D 10» This command allocates storage for 10 3D locations.

« **ADD3D** »

This command allows adding a 3D location to the 3D positions list.

Example: «ADD3D 10u 10u 0u» This command allocates storage for one location. Specify all axes locations using the same number format as MOVEX.

« **RUN3D** »

This command runs the stored sequence of 3D locations.

Example: «RUN3D»

« **DISIO** »

This command displays the setup of the TTL ports. DISIO must be used with a TTL port number (1-4).

Example :

« DISIO 1 »

Returns

«

status:

input

position:

Channel1

RisingEdge

»

« **CHAIO** »

This command allows setting the parameters of the 4 TTL IO ports.

Each TTL Port can be individually set to one of the follow states :

Disabled

Input Rising – responds to the rising edge of the incoming TTL signal

Input Falling – responds to the falling edge of the incoming TTL signal

Pulse Out Start – Provides a TTL pulse at the start of motion on the chosen axis

Pulse Out End – Provides a TTL pulse at the end of motion on the chosen axis

Pulse Out Top – Allows you to make the TTL signal at the best moment according to the shooting time you need.

Step Number Pulse – Provides a TTL pulse at the start of motion of the chosen axis on the specified step number

Gate pulse – Provides a TTL gate signal, switches on and off at the start of the specified steps

Output mode provides a TTL pulse at the start or end of motion of a chosen axis, or an output pulse on a specified step of motion on a given axis.

« CHAIO 1n » Disables TTL port 1

« CHAIO 1i1r » Sets TTL 1 as an input, triggering motion on the 1st axis on the rising edge

« CHAIO 1i2f » Sets TTL 1 as an input, triggering motion on the 2nd axis on the falling egde

« CHAIO 1o1s » Sets TTL 1 as an output, providing pulses at the start of motion of the first axis

« CHAIO 1o1t » Sets TTL 1 as an output, providing pulses at the best moment of motion of the first axis

« CHAIO 1o1e » Sets TTL 1 as an output, providing pulses at the end of motion of the first axis

« CHAIO 1o1n6 » Sets TTL 1 as an output, providing a pulse at the start of step number 6 motion on the first axis.

« CHAIO 1o1g5-10 » Sets TTL 1 as an output, provide a gate pulse controlled by the first axis, TTL output turns on at step 5 and turns off at step 10.

Format of disable commands :

« CHAIO [channel]d » where [channel] is TTL port 1-4

Format of Input commands :

« CHAIO [channel]i[axis][modifier] » where

[channel] is TTL port 1-4,

[axis] is the axis to trigger 1-3 and

[modifier] is either 'f' or 'r' designating triggering of the axis occurring on falling or rising edge of the TTL signal respectively.

Specific command have to be used for input TTL : the function to run a waveform with an input TTL is « RUTTL ». Use « CHAIA » for removing an input TTL CHAIO's command.

Format of the Output commands :

TTL pulses on every step:

« CHAIO [channel]o[axis][modifier] » where

[channel] is TTL port 1-4,

[axis] is the axis that generates TTL pulses and

[modifier] is either 's' or 'e' designating TTL pulses occurring at the start or end of motion respectively.

Or

TTL pulses on a given step number:

« CHAIO [channel]o[axis]n[step number] » where

[channel] is TTL port 1-4,

[axis] is the axis that generates TTL pulses and

'n' designating step number is to be provided

[step number] is an integer number within the step range as specified in ARBWF or SWF_n commands (eg for a waveform of 50 steps, the valid range of step numbers is 0-49)

Or

TTL Gate signal turned on and off at specified step numbers:

« CHAIO [channel]o[axis]g[start step number]-[end step number] » where

[channel] is TTL port 1-4,

[axis] is the axis that generates TTL pulses and

'g' designating TTL gate behaviour

[start step number] is an integer number within the step range as specified in ARBWF or SWF_n commands (eg for a waveform of 50 steps, the valid range of step numbers is 0-49)